

ADVISORY NOTICE

25_2025_XL 200 New software update

Dear customer,

We hereby inform you that the XL 200 device software has been updated to SW version 2025.01.

1. Purpose:

Application software ver.2025.01 released to include new features and to address issues reported in previous application software ver.2024.04A.

Embedded software ver. XL 200_2025_RS0.01 is release to address issues reported in previous embedded software ver. XL200_2024_RS0.13.

2. List of changes in Application Software ver. 2025.01:

New feature: Re-introduce option Use Test-wise Calibration in Calibration Policy Settings.

New Feature: Auto Span activity time is increased by ESW, Display Auto Span is in Progress.

New feature: When multiple barcoded bottles of same expiry are on-board; system will use the bottles as per FIFO.

SW version **2025.01** is already included in devices produced from SN 25P0384 onwards.

For SW versions, please contact your local service engineer or distributor.

The latest version of the SW is available on the [ERBA Knowledge Base](#) website.

Applied for:

INSTRUMENTS:

50006794	INS00001	XL-200 (with ISE)
50006770	INS00002	XL-200 (without ISE)

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno, 2nd June 2025

Ref. No.AN_25_2025