

ADVISORY NOTICE

27_2025_XL 200 New embedded software

Dear customer,

We hereby inform you that the XL 200 device Embedded software has been updated to SW version XL200_2024_RS0.13.

1. Purpose:

Embedded software ver. XL200_2024_RS0.13 is release to include improvement in previous embedded software ver. XL200_2023_RS0.09.

2. List of changes in Embedded software ver XL200_2024_RS0.13:

Improvement: Autospan Routine has been modified-(Performs autospan on Cuvette 1,15 & 29 and minimum gain value will be used for further reference).

Embedded SW version **XL200_2024_RS0.13** is already included in devices produced from SN 25P0340, 25P0342 onwards.

For Embedded SW version, please contact your local service engineer or distributor.
The latest version of the Embedded SW is available on the [ERBA Knowledge Base](#) website.

Applied for:

INSTRUMENTS:

50006794	INS00001	XL-200 (with ISE)
50006770	INS00002	XL-200 (without ISE)

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.




Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno, 2nd June 2025

Ref. No.AN_27_2025