

## ADVISORY NOTICE

28\_2025\_XL 640 New software update

Dear customer,

We hereby inform you that the XL 640 device software has been updated to SW version 2025.03A.

### 1. Purpose:

Application software 2025.03A is released to include new feature, improvements in previous Application software ver. 2024.05.

XL640+ Embedded software ESW2025\_RE0\_04 & FSW2025\_E0\_03 is released to include new feature and improvements in previous embedded software ESW2024\_RE0\_02 & FSW2024\_E0\_06.

### 2. List of changes in Application Software ver. 2025.03A:

Patient Reports - Translation Change: Caption is changed to Area, as suggested by Russian Team.

Re-introduce option Use Test-wise Calibration in Calibration Policy Settings.

New feature: When multiple barcoded bottles of same expiry are on-board; system will use the bottles as per FIFO

Facility to view Reagent Bottle History for particular Test/All Test.

SW version **2025.03A** is already included in devices produced from SN 645P103 onwards.

For SW versions, please contact your local service engineer or distributor.

The latest version of the SW is available on the [ERBA Knowledge Base](#) website.

Applied for:

### INSTRUMENTS:

50006795	INS00007	XL-640 (with ISE)
50006771	INS00008	XL-640 (without ISE)

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.



Brno, 2<sup>nd</sup> June 2025

Ref. No.AN\_28\_2025

