

ADVISORY NOTICE

34_2025_Erba Clean I recommended for EC 90

Dear customers,


We hereby inform you that the use of **Erba Clean I** (cat. No. EHL00042) is recommended for the internal decontamination process of the instrument EC 90.

The process of decontamination can be performed using either maximum 1% NaClO solution (as described in the User Manual) or the recommended Erba Clean I.

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.



Brno, 10th June 2025

Ref. No.AN_34_2025