

ADVISORY NOTICE

32_2025_XL 200 New software update

Dear customer,

We hereby inform you that the XL 200 device software has been updated to SW version 2025.01B.

1. Purpose:

Application software ver.2025.01B is released to include new feature in Previous ASW ver. 2025.01.

2. List of changes in Application Software ver. 2025.01B:

- New Feature: Provision for RUSSIA: Erba-RUS Product Login is added in software to display Logo Registered in Russia. (Refer Product Login document).
 Patient Reports - Translation Change: Caption is changed to Area.
 New Feature: CRR (Clinical Reportable Range) Scheme Revised:
 i. Pro-rate Tec-H Limit with (Decrease volume/ratio) for displaying Indicative values.
 ii. Do not Pro-rate Tec-L Limit with (Increase volume/ratio) for displaying Indicative values.

SW version **2025.01B** is already included in devices produced from SN 25P0393 onwards.

For SW versions, please contact your local service engineer or distributor.
 The latest version of the SW is available on the [ERBA Knowledge Base](#) website.

Applied for:


INSTRUMENTS:

50006794	INS00001	XL-200 (with ISE)
50006770	INS00002	XL-200 (without ISE)

In case you have any additional questions, please do not hesitate to contact us.

Best regards,
 Erba Lachema s.r.o.

Erba Lachema s.r.o.
 Karásek 2219/1d, 621 00 Brno
 DIČ: CZ26918846 -7-



Brno, 3th July 2025

Ref. No.AN_32_2025

