

## ADVISORY NOTICE

41\_2025\_ELite H5 CON

PRODUCT:	ELite H5 CON L	Cat. No. HEM00024
	ELite H5 CON N	Cat. No. HEM00025
	ELite H5 CON H	Cat. No. HEM00026
	ELite H5 CON	Cat. No. HEM00071

Dear customers,

We hereby inform you about a change in assayed values for **ELite H5 CON L/N/H, lot DH2507 / EH2507** due to inconsistency of QC results. This change is related to **H560** analysers with **“A” software versions only** and the modified values are highlighted in **bold**.

Please find the **new value sheet** and **inserts** on the website as usual. To verify you have the correct version of the value sheet, please check the *Date of Issue* which should be **24.07.2025**. To verify the insert, compare the insert values loaded in the instrument to the correct value sheet.

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

Erba Lachema s.r.o.  
Karásek 2219/1d, 621 00 Brno  
DIČ: CZ26918846 -7-



Brno, 24<sup>th</sup> July 2025

Ref. No.AN\_41\_2025

