

ADVISORY NOTICE

47_2025_New sw version for LAURA XL – 1.2.7.0

Dear customer,

We hereby inform you that there was a new software version released for LAURA XL automatic urine analyser. This sw version is 1.2.7.0 and is available from 10th September 2025. This sw update includes enhancements to the sediment evaluation in the presence of bubbles.

The new update includes:

- Improvement of bubble detection
- Best image selection for sediment evaluation

This new version is applied to new instruments production as of the release date. The upgrade is not mandatory for the installed instruments. It can however be beneficial when experiencing bubble presence.

Applied for:

INSTRUMENT: LAURA XL Cat. No:INS00065

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-



Brno, 10th September 2025

Ref. No.AN_47_2025