

ADVISORY NOTICE

56_2025_ Extension of declared values for glucose in Urinorm P

Dear customer,

We hereby inform you that the declared values for glucose (GLU) in Urinorm P have been extended to include the 100 mg/dL result. Based on our internal testing, the newly declared range for GLU in Urinorm is 100 – 1000 mg/dL / 5.5 – 55 mmol/L / 2+ – 4+. The instructions for use (IFU) for Urinorm are being updated accordingly. Please refer to this notice for the declared values until the revised IFU are available.

Applied for:

REAGENT:

URINORM

Cat. No. REG00053

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-



Brno, 21st October 2025

Ref. No.AN_56_2025