

## ADVISORY NOTICE

62\_2025\_CHEM-7\_New embedded software version v2024.C0.08

Dear customer,

We hereby inform you that the CHEM-7 device embedded software has been updated to SW version 2024.C0.08.

**1. Purpose:**

Embedded software ver. 2024.C0.08 for CHEM7 is released to address issue reported in previous embedded software ver.2024.C0.07.

**2. List of changes in Embedded Software ver. 2024.C0.08:**

Bug fixed: In PT Test Run, INR button & related action is not available

Embedded SW version **2024.C0.08** is already included in devices produced from SN S4050 onwards.

For SW versions, please contact your local service engineer or distributor.

The latest version of the SW is available on the [ERBA Knowledge Base](#) website. dd

Applied for:

**INSTRUMENTS:**

52000103	INS00014	CHEM-7 - ECZ
50005172	INS00014	CHEM-7
50008058	INS00014	CHEM-7 WITH INCUBATOR

In case you have any additional questions, please do not hesitate to contact us.

Best regards,  
Erba Lachema s.r.o.

**Erba Lachema s.r.o.**  
Karásek 2219/1d, 621 00 Brno  
DIČ: CZ26918846 -7-



Brno, 24<sup>th</sup> November 2025

Ref. No.AN\_62\_2025