

URINORM™

Control Urines



Urine analysis under control

Control Urines URINORM™



ANALYTES

- Specific gravity
- Leucocytes
- Nitrites
- Proteins
- pH
- Glucose
- Ketones
- Bilirubin
- Urobilinogen
- Blood

PRODUCT SPECIFICATION

- 2 level system of urine analysis control (normal and pathological levels)
- Packages: 2 x 3 vials of 15 ml
- Store at 2 - 8° C
- Stability of control urines following 1st opening: 1 month and/or 10 measurements
- Ready made solution for instant use
- 10 measurements with 1 vial
- Control urines to be used for PHAN® diagnostic strips and all LAURA® readers
- Quick results (within 1 min.)

Cat. nr.	Product	Qty
REG00053	Urinorm	2x3 bottles



Specifications are subject to change without prior notice.
Pictures used only for representation purposes - subject to change without prior notice.

Erba Lachema s.r.o.
 Karásek 2219/1d, 621 00 Brno, Czech Republic
 Tel.: +420 517 077 111
 E-mail: diagnostics@erbamannheim.com
www.erbalachema.com



Contact info / map

URINORM20180901EN