

# MAGO 4

## Infectious diseases



- Wide range of ELISA tests for infectious disease diagnostics
- Similar test protocols for individual product groups
- MAGO 4 pre-programmed test procedures
- Ready-to-use reagents, common product group components
- Capture methods for IgM antibody detection

# Infectious diseases

Code	Product	Packaging
967-821	CMV Screen	96 tests
967-820	CMV IgG	96 tests
967-830	CMV IgM	96 tests
967-837	CMV IgG Avidity	48 tests
720-600	Is EBV-VCA IgG	96 tests
720-610	Is EBV-VCA IgM	96 tests
720-640	Is EBV-EA-D IgG	96 tests
720-650	Is EBV-EA-D IgM	96 tests
720-620	Is EBNA-1 IgG	96 tests
720-630	Is EBNA-1 IgM	96 tests
967-840	Herpes 1 IgG	96 tests
967-850	Herpes 2 IgG	96 tests
967-845	Herpes 1+2 IgG	96 tests
967-835/M1	Herpes 1 IgM	96 tests
967-835/M2	Herpes 2 IgM	96 tests
967-835	Herpes 1+2 IgM	96 tests
967-800	Toxoplasma IgG	96 tests
967-810	Toxoplasma IgM	96 tests
967-805	Toxoplasma IgA	96 tests
967-807	Toxoplasma IgG Avidity	48 tests
961-100/R	Syphilis Screen Recombinant	96 tests
961-150	Treponema IgG	96 tests
961-151	Treponema IgM	96 tests
<b>Supplementary reagents</b>		
DBS-300	Sorbent (RF neutralisation)	120 tests

Complete list of products in Delta Biologicals catalogue.

- Automated processing of ELISA, IFA and agglutination tests in one batch
- 104 samples, barcodes, disposable tips
- Up to 12 ELISA or 16 IFA simultaneously in one batch
- Levey-Jennings charts and Westgard rules on board

