

ADVISORY NOTICE

13_2026_XL-1000_New software update

Dear customer,

We hereby inform you that the XL-1000 device application software has been updated to ASW version 2025.03B.

1. Purpose:

Application software ver.2025.03B released to include new feature in previous application software ver.2025.03A

2. List of changes in Application Software ver. 2025.03B:

New feature: Database performance was enhanced by relocating unused sample records to the archive.

Application SW version **2025.03B** is already included in devices produced from SN 115P015 onwards.

For SW versions, please contact your local service engineer or distributor.
The latest version of the SW is available on the [ERBA Knowledge Base](#) website.

Applied for:

INSTRUMENTS:

50006796	INS00009	XL-1000 (with ISE) plus
50006985	INS00010	XL-1000 (without ISE) plus

In case you have any additional questions, please do not hesitate to contact us.

Best regards,
Erba Lachema s.r.o.

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-



Brno, 12th March 2026

Ref. No.AN_13_2026