

ADVISORY NOTICE

24_2026_Total bilirubin in neonatal samples

Product: Bilirubin total DCA, **BIL T DCA 330**
Ref Code: XSYS0087

Dear customers,

With a mission to continuously enhance user experience and satisfaction, we would like to inform you about release of new test settings for BIL T DCA reagent (ref XSYS0087), intended for use for bilirubin total measurements in neonatal samples.

Recent studies demonstrated significant improvement in reagent performance and received patient results across neonatal samples while using these test settings.

To be able to use this kit for neonatal samples on Erba XL analyzers please upload test setting file named BTNEO, available in default parameters settings and Carry over pairs version 57, which can be downloaded from Erba Knowledge base.

In case of any questions, please do not hesitate to contact your local ERBA team.

Olga Kazakova

Best regards,

Olga Kazakova

Global Application Support Manager



Mailto: o.kazakova@erba.com

Brno, 22nd July 2026

Ref. No.AN_24_2026



Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, Czech Republic
Identification number: 269 18 846. Tax identification number: CZ26918846
Incorporated in the Commercial Register maintained by the Regional Court in Brno, Section C, insert 45458
Tel.: +420 517 077 111, e-mail: diagnostics@erba.com, www.erba.com