

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

Substance / mixture	Erba Chrom Protein C_R3
Number	mixture
UFI	EHL00009
Other mixture names	Y9AV-XW3W-XJ56-Y3KQ
R3: 1 x 5 ml PC Diluent	

1.2. Relevant identified uses of the substance or mixture and uses advised against**Mixture's intended use**

Reagent 3 of the kit Erba Chrom Protein C assay is intended for the quantitative determination of protein C in human plasma using a chromogenic assay method.

Main intended use

PC-MED-OTH Other medical devices

Secondary uses

PC-TEC-19 Reagents and laboratory chemicals

The use descriptors

PC 21 Laboratory chemicals

Mixture uses advised against

not available

1.3. Details of the supplier of the safety data sheet**Manufacturer**

Name or trade name	Erba Lachema s.r.o.
Address	Karásek 2219/1d , Brno, 62100
	Czech Republic
Identification number (CRN)	26918846
VAT number	CZ26918846
Phone	+420 517 077 111
Email	msds@erba.com
Web address	www.erba.com

Competent person responsible for the safety data sheet

Name	Erba Lachema s.r.o.
Email	msds@erba.com

1.4. Emergency telephone number

European emergency number: 112 112

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****Classification of the mixture in accordance with Regulation (EC) No 1272/2008**

The mixture is classified as dangerous.

Repr. 1B, H360D

Most serious adverse effects on human health and the environment

May damage the unborn child.

2.2. Label elements**Hazard pictogram****Signal word**

Danger

Hazardous substances

imidazole

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

Hazard statements

H360D May damage the unborn child.

Precautionary statements

P308+P313 IF exposed or concerned: Get medical advice/attention.

Supplemental information

EUH032 Contact with acids liberates very toxic gas.

2.3. Other hazards

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605. Mixture does not contain any substance meet the criteria for PBT or vPvB in accordance with Annex XIII of Regulation (EC) No. 1907/2006 (REACH) as amended. Does not contain any PMT or vPvM components.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Mixture of substances and additives specified below.

Mixture contains these hazardous substances and substances with the highest permissible concentration in the working environment

Identification numbers	Substance name	Content in % weight	Classification according to Regulation (EC) No 1272/2008	Note
Index: 011-004-00-7 CAS: 26628-22-8 EC: 247-852-1	sodium azide	<1	Acute Tox. 2, H300 Aquatic Acute 1, H400 Aquatic Chronic 1, H410 EUH032	1
Index: 613-319-00-0 CAS: 288-32-4 EC: 206-019-2 Registration number: 01-2119485825-24-x	imidazole	<1	Acute Tox. 4, H302 Skin Corr. 1C, H314 Repr. 1B, H360D	

Notes

1 A substance for which exposure limits are set.

Full text of all classifications and hazard statements is given in the section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures**

Take care of your own safety. If any health problems are manifested or if in doubt, inform a doctor and show him information from this safety data sheet.

If inhaled

Discontinue the exposure, remove casualty to fresh air, keep at rest and seek medical advice.

If on skin

Remove contaminated clothes. Wash the affected area with plenty of water, lukewarm if possible. Soap, soap solution or shampoo should be used if there is no skin injury.

If in eyes

Rinse an open eye (hold eyelids with fingers) with plenty of water for about 15 minutes, transfer casualty to a specialist.

If swallowed

Rinse mouth with water, drink 1/2 l of lukewarm water, DO induce vomiting. Seek medical attention immediately.

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

- 4.2. Most important symptoms and effects, both acute and delayed**
If inhaled
Not expected.
If on skin
Irritation, redness.
If in eyes
Not expected.
If swallowed
Irritation, nausea, stomach pain, vomiting, tightness in the chest.
- 4.3. Indication of any immediate medical attention and special treatment needed**
Symptomatic treatment.

SECTION 5: Firefighting measures

- 5.1. Extinguishing media**
Suitable extinguishing media
The mixture is not flammable. The measures should be adapted to burning substances in the surrounding area.
Unsuitable extinguishing media
For mixture no unsuitable extinguishing media are known.
- 5.2. Special hazards arising from the substance or mixture**
In the event of fire, carbon monoxide, carbon dioxide and other toxic gases may arise. Inhalation of hazardous degradation (pyrolysis) products may cause serious health damage.
- 5.3. Advice for firefighters**
Self-Contained Breathing Apparatus (SCBA) with a chemical protection suit only where personal (close) contact is likely. Use a self-contained breathing apparatus and full-body protective clothing. Do not allow run-off of contaminated fire extinguishing material to enter drains or surface and ground water.

SECTION 6: Accidental release measures

- 6.1. Personal precautions, protective equipment and emergency procedures**
Pregnant women should strictly avoid inhalation and contact with the product. Use personal protective equipment, see Section 8. Observe the principles of work safety in chemical laboratories. Do not eat, drink or smoke.
- 6.2. Environmental precautions**
Do not discharge into the drains, surface waters and groundwater.
- 6.3. Methods and material for containment and cleaning up**
Absorb spilled agent with a suitable inert material (sand, earth, vapex) and store contaminated material in containers for collection of hazardous waste. For waste disposal, see Section 13.
- 6.4. Reference to other sections**
See the Section 7, 8 and 13.

SECTION 7: Handling and storage

- 7.1. Precautions for safe handling**
Avoid contact during pregnancy and while nursing. Observe the normal operating procedures for handling chemical substances and mixtures. Use personal protective equipment, see Section 8.
- 7.2. Conditions for safe storage, including any incompatibilities**
Store in tightly closed containers in cold, dry and well ventilated areas designated for this purpose.
Storage temperature min 2 °C, max 8 °C
- 7.3. Specific end use(s)**
The reagent of the kit is designed for in vitro diagnostic use.

SECTION 8: Exposure controls/personal protection

- 8.1. Control parameters**
The mixture contains substances for which occupational exposure limits are set.

Erba Chrom Protein C_R3

Creation date 20th February 2015 Version 4.0
Revision date 16th September 2026

Czech Republic

Government Regulation 473/2025 Coll.

Substance name (component)	Type	Value
sodium azide (CAS: 26628-22-8)	PEL	0,1 mg/m ³
	NPK-P	0,3 mg/m ³

Notes

Skin penetration is significantly involved during exposure.

Irritating to mucous membranes (eyes, respiratory system) and skin.

European Union

Commission Directive 2000/39/EC

Substance name (component)	Type	Value
sodium azide (CAS: 26628-22-8)	OEL 8 hours	0,1 mg/m ³
	OEL 15 minutes	0,3 mg/m ³

Notes

Skin.

United Kingdom

EH40/2005 Workplace exposure limits (Fourth Edition 2020)

Substance name (component)	Type	Value
sodium azide (as NaN ₃) (CAS: 26628-22-8)	WEL 8h	0,1 mg/m ³
	WEL 15min	0,3 mg/m ³

Notes

Can be absorbed through the skin. The assigned substances are those for which there are concerns that dermal absorption will lead to systemic toxicity.

8.2. Exposure controls

Follow the usual measures intended for health protection at work and especially for good ventilation. This can be achieved only by local suction or efficient general ventilation. Do not eat, drink and smoke during work. Wash your hands thoroughly with water and soap after work and before breaks for a meal and rest. Pregnant women should strictly avoid inhalation or skin contact.

Eye/face protection

Protective goggles.

Skin protection

Hand protection: Protective gloves resistant to the product. Contaminated skin should be washed thoroughly. Other protection: protective clothing.

Respiratory protection

If all workplace limits are observed and good ventilation is ensured, no special precautions necessary.

Thermal hazard

Not available.

Environmental exposure controls

Observe usual measures for protection of the environment, see Section 6.2.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	liquid
Colour	colourless
Odour	odourless
Melting point/freezing point	data not available
Boiling point or initial boiling point and boiling range	data not available
Flammability	data not available
Lower and upper explosion limit	data not available
Flash point	data not available
Auto-ignition temperature	data not available
Decomposition temperature	data not available
pH	data not available
Kinematic viscosity	data not available

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

Solubility in water	data not available
Solubility in fats	data not available
Partition coefficient n-octanol/water (log value)	data not available
Vapour pressure	data not available
Density and/or relative density	data not available
Relative vapour density	data not available
Particle characteristics	data not available

9.2. Other information

Evaporation rate	data not available
none	

SECTION 10: Stability and reactivity**10.1. Reactivity**

Under normal conditions of use and storage the mixture is stable.

10.2. Chemical stability

The product is stable under normal conditions.

10.3. Possibility of hazardous reactions

Contact with acids may liberate toxic gas.

10.4. Conditions to avoid

The product is stable and no degradation occurs under normal use.

10.5. Incompatible materials

Protect against strong acids, strong oxidizing agents.

10.6. Hazardous decomposition products

Toxic fumes.

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008**

Inhalation of solvent vapors above values exceeding exposure limits for working environment may result in acute inhalation poisoning, depending on the level of concentration and exposure time. No toxicological data is available for the mixture.

Acute toxicity

Based on the available data, the criteria for classification of the mixture are not met.

Erba Chrom Protein C_R3						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	ATE	11705 mg/kg				Calculation of value

imidazole						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	880 mg/kg		Mouse		
	LD ₅₀	970 mg/kg		Rat (Rattus norvegicus)		

sodium azide						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	27 mg/kg bw				
Inhalation	LC ₅₀	54 mg/m ³	4 hours	Rat		

Skin corrosion/irritation

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

Serious eye damage/irritation

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Respiratory or skin sensitisation

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Germ cell mutagenicity

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Carcinogenicity

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Reproductive toxicity

May damage the unborn child. Data for the components of the mixture are not available.

Toxicity for specific target organ - single exposure

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Toxicity for specific target organ - repeated exposure

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Aspiration hazard

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

11.2. Information on other hazards**Endocrine disrupting properties**

Based on the available data, the criteria for classification of the mixture are not met. Does not contain any components that may cause endocrine disruption for humans.

Other information

not available

SECTION 12: Ecological information**12.1. Toxicity**

Data for the mixture are not available. Based on the available data, the criteria for classification of the mixture are not met.

Acute toxicity

imidazole				
Parameter	Value	Exposure time	Species	Environment
EC ₅₀	231 mg/l	30 minutes		
EC ₅₀	341.5 mg/l	48 hours		

sodium azide				
Parameter	Value	Exposure time	Species	Environment
LC ₅₀	680 µg/l		Fish	Fresh water
EC ₅₀ /LC ₅₀	400 µg/l		Invertebrates	Fresh water

Erba Chrom Protein C_R3

Creation date 20th February 2015 Version 4.0
Revision date 16th September 2026

sodium azide				
Parameter	Value	Exposure time	Species	Environment
EC ₅₀ /LC ₅₀	150 µg/l		Invertebrates	Salt water
EC ₅₀ /LC ₅₀	348 µg/l		Algae	Fresh water
EC ₅₀ /LC ₅₀	5.6 mg/l		Microorganisms	
NOEC	30 µg/l		Microorganisms	

12.2. Persistence and degradability

No data are available for either the mixture or the components.

12.3. Bioaccumulative potential

No data are available for either the mixture or the components.

12.4. Mobility in soil

No data are available for either the mixture or the components.

12.5. Results of PBT and vPvB assessment

Based on the available data, the criteria for classification of the mixture are not met. Does not contain any PBT or vPvB components. Product does not contain any substance meeting the criteria for PBT or vPvB in accordance with the Annex XIII of Regulation (EC) No 1907/2006 (REACH) as amended.

12.6. Endocrine disrupting properties

Based on the available data, the criteria for classification of the mixture are not met. Does not contain any components that may cause endocrine disruption in the environment.

12.7. Other adverse effects

Negligible ecotoxicity.

SECTION 13: Disposal considerations**13.1. Waste treatment methods**

Hazard of environmental contamination; dispose of the waste in accordance with the local and/or national regulations. Any unused product and contaminated packaging should be put in labelled containers for waste collection and submitted for disposal to a person authorised for waste removal (a specialized company) that is entitled for such activity. Do not empty unused product in drainage systems. The product must not be disposed of with municipal waste. Perfectly cleaned containers can be submitted for recycling.

Waste management legislation

Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste, as amended. Decision 2000/532/EC establishing a list of wastes, as amended.

Waste type code

16 05 06* laboratory chemicals, consisting of or containing hazardous substances, including mixtures of laboratory chemicals

Packaging waste type code

15 01 10* packaging containing residues of or contaminated by hazardous substances

(*) - Hazardous waste according to Directive 2008/98/EC on hazardous waste

SECTION 14: Transport information**14.1. UN number or ID number**

not subject to transport regulations

14.2. UN proper shipping name

not relevant

14.3. Transport hazard class(es)

not relevant

14.4. Packing group

not relevant

14.5. Environmental hazards

not relevant

14.6. Special precautions for user

Reference in the Sections 4 to 8.

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

14.7. Maritime transport in bulk according to IMO instruments

Not transported.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18th December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing the European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No. 793/93 and Commission Regulation (EC) No. 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, as amended. REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16th December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No. 1907/2006, as amended. Commission Regulation (EU) 2020/878 of 18 June 2020 amending Annex II to Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).

15.2. Chemical safety assessment

Assessment was not carried out.

SECTION 16: Other information**A list of standard risk phrases used in the safety data sheet**

EUH032	Contact with acids liberates very toxic gas.
H300	Fatal if swallowed.
H302	Harmful if swallowed.
H314	Causes severe skin burns and eye damage.
H360D	May damage the unborn child.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.

Guidelines for safe handling used in the safety data sheet

P308+P313 IF exposed or concerned: Get medical advice/attention.

Other important information about human health protection

The product must not be - unless specifically approved by the manufacturer/importer - used for purposes other than as per the Section 1. The user is responsible for adherence to all related health protection regulations.

Key to abbreviations and acronyms used in the safety data sheet

Acute Tox.	Acute toxicity
ADR	Agreement concerning the international carriage of dangerous goods by road
Aquatic Acute	Hazardous to the aquatic environment
Aquatic Chronic	Hazardous to the aquatic environment (chronic)
ATE	Acute toxicity estimate
BCF	Bioconcentration Factor
CAS	Chemical Abstracts Service
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substance and mixtures
EC	Identification code for each substance listed in EINECS
EC ₅₀	Concentration of a substance when it is affected 50 % of the population
EINECS	European Inventory of Existing Commercial Chemical Substances
EmS	Emergency Response Procedures for Ships Carrying Dangerous Goods
EU	European Union
EuPCS	European Product Categorisation System
IATA	International Air Transport Association
IBC	International Code For The Construction And Equipment of Ships Carrying Dangerous Chemicals
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization

Erba Chrom Protein C_R3

Creation date	20th February 2015	Version	4.0
Revision date	16th September 2026		

INCI	International Nomenclature of Cosmetic Ingredients
ISO	International Organization for Standardization
IUPAC	International Union of Pure and Applied Chemistry
LC ₅₀	Lethal concentration of a substance in which it can be expected death of 50% of the population
LD ₅₀	Lethal dose of a substance in which it can be expected death of 50% of the population
log K _{ow}	Octanol-water partition coefficient
NOEC	No observed effect concentration
NPK	Maximum admissible concentration
OEL	Occupational Exposure Limits
PBT	Persistent, bioaccumulative and toxic
PEL	Permissible Exposure Limit
PMT	Persistent, mobile and toxic
ppm	Parts per million
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
Repr.	Reproductive toxicity
RID	Regulation concerning the International Carriage of Dangerous Goods by Rail
Skin Corr.	Skin corrosion
UN number	Four-figure identification number of the substance or article taken from the UN Model Regulations
UVCB	Substances of unknown or variable composition, complex reaction products or biological materials
VOC	Volatile organic compounds
vPvB	Very persistent and very bioaccumulative
vPvM	Very persistent and very mobile

Training guidelines

Inform the personnel about the recommended ways of use, mandatory protective equipment, first aid and prohibited ways of handling the product.

Recommended restrictions of use

This compound is designed for professional use. It should not be used for purposes other than those described in Section 1.2.

Information about data sources used to compile the Safety Data Sheet

REGULATION (EC) No. 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL (REACH) as amended.
REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Data from the manufacturer of the substance / mixture, if available - information from registration dossiers.

The changes (which information has been added, deleted or modified)

The version 4.0 replaces the SDS version from Thursday, 19 October 2023. Changes were made in sections 2, 11, 12, 13 and 16.

Statement

The safety data sheet provides information aimed at ensuring safety and health protection at work and environmental protection. The provided information corresponds to the current status of knowledge and experience and complies with valid legal regulations. The information should not be understood as guaranteeing the suitability and usability of the product for a particular application.