

# SYSPACKS™

Reagents for  
ERBA XL Analysers



 Erba Lachema s.r.o.  
Karásek 1d, 621 00 Brno, Czech Republic  
Tel.: +420 517 077 111, fax: +420 517 077 077  
E-mail: [diagnostics@erbalachema.com](mailto:diagnostics@erbalachema.com)  
[www.erbalachema.com](http://www.erbalachema.com)



Contact info / map

Ready-to-Use  
Reliability

## CLINICAL CHEMISTRY

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0001	ALBUMIN	ALB 440	R1:10 x 44 ml			BCG
XSYS0002	ALKALINE PHOSPHATASE	ALP 110	R1: 2 x 44 ml	R2: 2 x 11 ml		optim. IFCC (AMP buffer)
XSYS0003	ALPHA AMYLASE (SINGLE REAGENT)	AMY 110	R1: 5 x 22 ml			CNP-G, single reagent
XSYS0017	ALT/GPT	ALT/GPT 330	R1: 6 x 44 ml	R2: 3 x 22 ml		mod. IFCC
XSYS0074	ALT/GPT	ALT/GPT 564 XL-1000	R1: 6 x 72 ml	R2: 6 x 22 ml		mod. IFCC
XSYS0016	AST/GOT	AST/GOT 330	R1: 6 x 44 ml	R2: 3 x 22 ml		mod. IFCC
XSYS0073	AST/GOT	AST/GOT 564 XL-1000	R1: 6 x 72 ml	R2: 6 x 22 ml		mod. IFCC
XSYS0028	BILIRUBIN DIRECT	BIL D 330	R1: 6 x 44 ml	R2: 3 x 22 ml		Walter-Gerarde
XSYS0023	BILIRUBIN TOTAL	BIL T 330	R1: 6 x 44 ml	R2: 3 x 22 ml		Walter-Gerarde
XSYS0086	BILIRUBIN DIRECT DCA	BIL D DCA 330	R1: 6 x 44 ml	R2: 3 x 22 ml		with DCA
XSYS0087	BILIRUBIN TOTAL DCA	BIL T DCA 330	R1: 6 x 44 ml	R2: 3 x 22 ml		with DCA
XSYS0007	CALCIUM	CA 120	R1: 10 x 12 ml			Arsenaso III
XSYS0022	CREATINE KINASE	CK 110	R1: 2 x 44 ml	R2: 2 x 11 ml		optim. IFCC, DGKCH
XSYS0029	CREATINE KINASE MB	CK MB 110	R1: 2 x 44 ml	R2: 2 x 11 ml		immuno-inhibition
XSYS0024	CREATININE	CREA 275	R1: 5 x 44 ml	R2: 5 x 11 ml		Jaffé
XSYS0076	CREATININE	CREA 564 XL-1000	R1: 6 x 72 ml	R2: 6 x 22 ml		Jaffé
XSYS0085	CREATININE ENZYMATIC	CREA ENZ 200	R1: 5 x 30 ml	R2: 5 x 10 ml		enzymatic
XSYS0011	GAMMAGLUTAMYLTRANSFERASE	GGT 110	R1: 2 x 44 ml	R2: 2 x 11 ml		Persijn van der Slik
XSYS0077	GAMMAGLUTAMYLTRANSFERASE	GGT 330 XL-1000	R1: 6 x 44 ml	R2: 3 x 22 ml		Persijn van der Slik
XSYS0012	GLUCOSE	GLU 440	R1:10 x 44 ml			GOD,POD
XSYS0069	GLUCOSE	GLU 576 XL-1000	R1: 8 x 72 ml			GOD,POD
XSYS0043	HDL DIRECT	HDL C 160	R1: 4 x 30 ml	R2: 4 x 10 ml		direct method with PVS,PEGME
XSYS0078	HDL DIRECT	HDL C 360 XL-1000	R1: 4 x 68 ml	R2: 4 x 22 ml		direct method with PVS,PEGME
XSYS0008	CHLORIDE	CL 120	R1: 10 x 12 ml			thiocyanate
XSYS0009	CHOLESTEROL	CHOL 440	R1:10 x 44 ml			CHOD,PAP
XSYS0070	CHOLESTEROL	CHOL 576 XL-1000	R1: 8 x 72 ml			CHOD,PAP
XSYS0049	IRON	FE 125	R1: 4 x 25 ml	R2: 2 x 12,5 ml	Standard: 2 x 2 ml	ferrozine
XSYS0013	LACTATEDEHYDROGENASE-P	LDH 110	R1: 2 x 44 ml	R2: 2 x 11 ml		optim. DGKCH (pyruvate)
XSYS0044	LDL DIRECT	LDL C 80	R1: 2 x 30 ml	R2: 2 x 10 ml		direct method with PVS,PEGME
XSYS0081	LIPASE	LIP 110	R1: 2 x 44 ml	R2: 2 x 11 ml		enzymatic, colorimetric
XSYS0040	MAGNESIUM	MG 88	R1: 2 x 44 ml			xylidyl blue
XSYS0027	MICROPROTEIN	MP 120	R1: 10 x 12 ml		Standard: 1x 5 ml	pyrogallol red
XSYS0015	PHOSPHORUS	PHOS 120	R1: 10 x 12 ml			UV molybdate
XSYS0018	TOTAL PROTEIN	TP 440	R1:10 x 44 ml			biuret, single reagent
XSYS0041	TRIGLYCERIDES	TG 440	R1:10 x 44 ml			GPO,PAP, single reagent
XSYS0071	TRIGLYCERIDES	TG 576 XL-1000	R1: 8 x 72 ml			GPO,PAP, single reagent
XSYS0050	UIBC	UIBC 125	R1: 4 x 25 ml	R2: 2 x 12,5 ml	Standard: 2 x 2 ml	ferrozine
XSYS0020	UREA	UREA 275	R1: 5 x 44 ml	R2: 5 x 11 ml		Urease, GLDH
XSYS0075	UREA	UREA 564 XL-1000	R1: 6 x 72 ml	R2: 6 x 22 ml		Urease, GLDH
XSYS0042	URIC ACID (SINGLE REAGENT)	UA 440	R1:10 x 44 ml			enzymatic, single reagent
XSYS0072	URIC ACID (SINGLE REAGENT)	UA 576 XL-1000	R1: 8 x 72 ml			enzymatic, single reagent

## IMMUNOASSAY

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0046	ANTISTREPTOLYSIN (O)	ASO	R1: 2 x 40 ml	R2: 2 x 10 ml		latex enhanced immunoturbidimetry
XSYS0047	C-REACTIVE PROTEIN	CRP	R1: 2 x 40 ml	R2: 2 x 10 ml		immunoturbidimetry
XSYS0084	C- REACTIVE PROTEIN HIGH SENSITIVE	CRP HS	R1: 2 x 40 ml	R2: 2 x 8 ml		immunoturbidimetry
XSYS0054	HbA1C	HbA1C	R1: 1 x 24 ml	R2a: 1 x 8 ml R2b: 1 x 4 ml	R3: 2 x 50 ml	latex enhanced immunoturbidimetry
XSYS0083	MICROALBUMIN	MAL	R1: 2 x 30 ml	R2: 1 x 10 ml		immunoturbidimetry
XSYS0048	RHEUMATOID FACTOR	RF	R1: 2 x 40 ml	R2: 2 x 8 ml		immunoturbidimetry

## CALIBRATION AND CONTROL MATERIALS

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0034	XL MULTICAL	XL MULTICAL	R1: 4x3 ml			certified human multicalibrator
XSYS0061	HDL/LDL CALIBRATOR	HDL/LDL CAL	R1: 2x1 ml			certified human calibrator for HDL and LDL cholesterol
BLT00080	ERBA NORM	ERBA NORM	R1: 4x5 ml	Water 1x20 ml		certified human multicontrol sera
BLT00081	ERBA PATH	ERBA PATH	R1: 4x5 ml	Water 1x20 ml		certified human multicontrol sera
XSYS0051	ASO CALIBRATOR	ASO CAL SH	R1: 1 x 1 ml			certified human calibrator for ASO
BLT20004	ASO CONTROL	ASO CON	R1: 1 x 1 ml			certified human control serum for ASO
XSYS0053	CRP CALIBRATOR	CRP CAL SH	R1: 1 x 1 ml			certified human calibrator for CRP
XSYS0088	CRP-HS CALIBRATOR	CRP-HS CAL	R1: 1 x 1 ml			certified human calibrator for CRP-HS
BLT20013	CRP CONTROL	CRP CON H	R1: 1 x 1 ml			certified human control serum for CRP
BLT20014	CRP CONTROL	CRP CON L	R1: 1 x 1 ml			certified human control serum for CRP
XSYS0057	HbA1c CALIBRATOR SET	HbA1c CAL SET	R1: 4 x 0.25 ml			certified human calibrator for HbA1C
XSYS0056	HbA1c CONTROL	HbA1c CON H	R1: 1 x 0.25 ml			certified human control serum for HbA1C
XSYS0055	HbA1c CONTROL	HbA1c CON L	R1: 1 x 0.25 ml			certified human control serum for HbA1C
XSYS0052	RF CALIBRATOR	RF CAL SH	R1: 1 x 1 ml			certified human calibrator for RF
BLT20039	RF CONTROL	RF CON	R1: 1 x 1 ml			certified human control serum for RF

## SOLUTIONS FOR ISE

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0064	XL CALIBRANT B	XL CAL A	R1: 3x500 ml			water standard
XSYS0065	XL CALIBRANT A	XL CAL B	R1: 1x52 ml			water standard
XSYS0062	XL ISE CLEANING	XL ISE CLEAN	R1: 1x52 ml			kinase based cleaner
XSYS0063	XL ISE 4 CLEANING	XL ISE 4 CLEAN	R1: 1x90 ml	R2: 6 x 0.3 g		pepsin based cleaner
XSYS0067	XL URINE DILUENT	XL URINE DIL	R1: 1x52 ml			urine diluent

## WASHING SOLUTIONS

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0066	XL WASH	XL WASH	R1: 4x100 ml			wash solution for daily washing
XSYS0058	XL AUTOWASH AC	XL AUTOWASH AC	R1: 1x500 ml			acid wash solution for occasional washing
XSYS0059	XL AUTOWASH AL	XL AUTOWASH AL	R1: 1x500 ml			acid wash solution for occasional washing
XSYS0082	XL AUTOWASH AC/AL	XL AUTOWASH AC/ AL	R1: 5x44 ml	R1:5x44 ml		acid and alkali solution for occasional washing